

For Research Use Only

# Flotillin 2 Polyclonal antibody

Catalog Number: 11446-1-AP

2 Publications



## Basic Information

Catalog Number:

11446-1-AP

Size:

150ul, Concentration: 287 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1996

GenBank Accession Number:

BC017292

GeneID (NCBI):

2319

UNIPROT ID:

Q14254

Full Name:

flotillin 2

Calculated MW:

379 aa, 42 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

Flotillins are membrane-associated proteins that contain two homologous isoforms, flotillin 1 (FLOT1) and flotillin 2 (FLOT2). They are highly conserved and widely expressed proteins that commonly used as markers of lipid rafts. Flotillins are involved in various signaling pathways, cell adhesion, membrane trafficking and axonal growth.

## Notable Publications

Author	Pubmed ID	Journal	Application
Min Cheng	39265655	Antiviral Res	
Oindrila Paul	38771750	PLoS One	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

