

For Research Use Only

Flotillin 2 Polyclonal antibody

Catalog Number: 11446-1-AP



Basic Information

Catalog Number: 11446-1-AP	GenBank Accession Number: BC017292	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 287 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2319	
Source: Rabbit	Full Name: flotillin 2	
Isotype: IgG	Calculated MW: 379 aa, 42 kDa	
Immunogen Catalog Number: AG1996		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Flotillins are membrane-associated proteins that contain two homologous isoforms, flotillin 1 (FLOT1) and flotillin 2 (FLOT2). They are highly conserved and widely expressed proteins that commonly used as markers of lipid rafts. Flotillins are involved in various signaling pathways, cell adhesion, membrane trafficking and axonal growth.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data