

For Research Use Only

RAX2 Polyclonal antibody

Catalog Number: 11466-1-AP **1 Publications**



Basic Information

Catalog Number: 11466-1-AP	GenBank Accession Number: BC018709	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1200 µg/ml by Nanodrop and 513 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 84839	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: retina and anterior neural fold homeobox 2	
Isotype: IgG	Calculated MW: 184 aa, 20 kDa	
Immunogen Catalog Number: AG2046	Observed MW: 22-27 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells,
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Rui Su	32015336	Cell Death Dis	WB

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 11466-1-AP (RAX2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.