

For Research Use Only

RAX2 Polyclonal antibody

Catalog Number: 11466-1-AP

1 Publications



Basic Information

Catalog Number:

11466-1-AP

Size:

150ul, Concentration: 1200 ug/ml by Nanodrop and 513 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2046

GenBank Accession Number:

BC018709

GeneID (NCBI):

84839

UNIPROT ID:

Q96IS3

Full Name:

retina and anterior neural fold homeobox 2

Calculated MW:

184 aa, 20 kDa

Observed MW:

22-27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : Jurkat cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Rui Su	32015336	Cell Death Dis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 11466-1-AP (RAX2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.