

For Research Use Only

ERAL1 Polyclonal antibody

Catalog Number: 11478-1-AP

Featured Product

12 Publications



Basic Information

Catalog Number:

11478-1-AP

Size:

150ul , Concentration: 347 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2019

GenBank Accession Number:

BC019094

GeneID (NCBI):

26284

UNIPROT ID:

O75616

Full Name:

Era G-protein-like 1 (E. coli)

Calculated MW:

437 aa, 48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF, ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Agata Rozanska	28512204	Biochem J	WB
Tian Luo	28553621	Front Cell Infect Microbiol	WB,IF
Di Hu	30877431	Acta Neuropathol	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

