

For Research Use Only

STAC Polyclonal antibody

Catalog Number: 11480-1-AP



Basic Information

Catalog Number: 11480-1-AP	GenBank Accession Number: BC020221	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop and 433 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 6769	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: SH3 and cysteine rich domain	
Isotype: IgG	Calculated MW: 402 aa, 45 kDa	
Immunogen Catalog Number: AG2059	Observed MW: 45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue, PC-3 cells, human liver tissue
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 11480-1-AP (STAC antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.