

For Research Use Only

MFSD10 Polyclonal antibody, PBS Only

Catalog Number: 11518-1-PBS



Basic Information

Catalog Number:

11518-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2080

GenBank Accession Number:

BC014979

GeneID (NCBI):

10227

UNIPROT ID:

Q14728

Full Name:

major facilitator superfamily domain containing 10

Calculated MW:

455 aa, 48 kDa

Observed MW:

48-52 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 11518-1-AP (MFSD10 Antibody) at dilution of 1:600 incubated at 4 degree celsius over night. This data was developed using the same antibody clone with 11518-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 11518-1-AP (MFSD10 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11518-1-PBS in a different storage buffer formulation.