

For Research Use Only

DNTTIP1 Polyclonal antibody, PBS Only

Catalog Number: 11637-1-PBS



Basic Information

Catalog Number:

11637-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2213

GenBank Accession Number:

BC024290

GeneID (NCBI):

116092

UNIPROT ID:

Q9H147

Full Name:

deoxynucleotidyltransferase, terminal, interacting protein 1

Calculated MW:

329 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, IP, Indirect ELISA

Species Specificity:

human

Background Information

DNTTIP1, also named as Deoxynucleotidyltransferase terminal-interacting protein 1, is a 329 amino acid protein, which contains a Nuclear localization signal and localizes in the nuclear. DNTTIP1 also acts as a transcriptional regulator, which binding to the consensus sequence 5'-GNTGCATG-3' following an AT-tract. Associates with RAB20 promoter and positively regulates its transcription.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cell were subjected to SDS PAGE followed by western blot with 11637-1-AP (DNTTIP1 Antibody) at dilution of 1:500. This data was developed using the same antibody clone with 11637-1-PBS in a different storage buffer formulation.



IP result of anti-DNTTIP1 (IP:11637-1-AP, 4ug; Detection:11637-1-AP 1:300) with HeLa cells lysate 1200ug. This data was developed using the same antibody clone with 11637-1-PBS in a different storage buffer formulation.

Immunofluorescent analysis of (4% PFA) fixed HeLa cells using DNTTIP1 antibody (11637-1-AP) at dilution of 1:400 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002), CL594-Phalloidin (red). This data was developed using the same antibody clone with 11637-1-PBS in a different storage buffer formulation.