

For Research Use Only

ETV3 Polyclonal antibody, PBS Only

Catalog Number:11721-1-PBS



Basic Information

Catalog Number:

11721-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2301

GenBank Accession Number:

BC022868

GeneID (NCBI):

2117

UNIPROT ID:

P41162

Full Name:

ets variant 3

Calculated MW:

512 aa, 57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

ETV3 is a transcriptional repressor that belongs to the ETS family of transcription factors, which bind a common DNA sequence core motif (GGAA) via their ETS domains. ETV3 is preferentially expressed in mature DC including tissue-derived migratory DC (migDC) and facilitates their homeostatic maturation and CCR7-dependent migration. Phosphorylation of ETV3 regulated transcription by preventing its binding to DNA at promoters for several thousand genes, including some involved in negative feedback regulation of itself and of upstream signals (PMID: 36543959, 41678619). Western and IP analysis detected the specific band with an apparent MW of 57 kDa.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

