

For Research Use Only

PPP1R1C Polyclonal antibody

Catalog Number: 11757-1-AP



Basic Information

Catalog Number:

11757-1-AP

Size:

150ul , Concentration: 213 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2353

GenBank Accession Number:

BC017943

GeneID (NCBI):

151242

UNIPROT ID:

Q8WV17

Full Name:

protein phosphatase 1, regulatory (inhibitor) subunit 1C

Calculated MW:

13 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

