

For Research Use Only

FCN1 Polyclonal antibody

Catalog Number: 11775-1-AP

3 Publications



Basic Information

Catalog Number:

11775-1-AP

Size:

150ul, Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2444

GenBank Accession Number:

BC020635

GeneID (NCBI):

2219

UNIPROT ID:

O00602

Full Name:

ficolin (collagen/fibrinogen domain containing) 1

Calculated MW:

326 aa, 35 kDa

Observed MW:

33-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB: human plasma, K-562 cells, human liver tissue, mouse liver tissue, rat liver tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Muhammad Khan	36389725	Front Immunol	IHC
Kohei Kusakari	36052069	Front Immunol	WB
Fangfang Zhong	39139822	Front Genet	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11775-1-AP (FCN1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.