

For Research Use Only

AFTPH Polyclonal antibody

Catalog Number: 11785-1-AP



Basic Information

Catalog Number:

11785-1-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2402

GenBank Accession Number:

BC022247

GeneID (NCBI):

54812

UNIPROT ID:

Q6ULP2

Full Name:

aftiphilin

Calculated MW:

937 aa, 102 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Aftiphilin is a component of the clathrin machinery in neurons, binding to of clathrin adaptor protein 1 and 2. Western assay showed that the protein emigrate as a 126 kDa protein in SDS PAGE, with larger molecular weight than the predicted size (99 and 102 kDa for the two isoforms).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

