

For Research Use Only

BUD31 Polyclonal antibody, PBS Only

Catalog Number: 11798-1-PBS

Featured Product



Basic Information

Catalog Number:

11798-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2373

GenBank Accession Number:

BC022821

GeneID (NCBI):

8896

UNIPROT ID:

P41223

Full Name:

BUD31 homolog (S. cerevisiae)

Calculated MW:

144 aa, 17 kDa

Observed MW:

17-20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 11798-1-AP (BUD31 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11798-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer using 11798-1-AP (BUD31 antibody) at dilution of 1:100 (under 10x lens). This data was developed using the same antibody clone with 11798-1-PBS in a different storage buffer formulation.



IP result of anti-BUD31 (IP:11798-1-AP, 4ug; Detection:11798-1-AP 1:1500) with mouse brain tissue lysate 3440ug. This data was developed using the same antibody clone with 11798-1-PBS in a different storage buffer formulation.