

For Research Use Only

RABGGTB Polyclonal antibody

Catalog Number: 11801-1-AP



Basic Information

Catalog Number: 11801-1-AP	GenBank Accession Number: BC020790	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 µg/ml by Nanodrop and 267 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 5876	Recommended Dilutions: WB 1:1000-1:3000
Source: Rabbit	Full Name: Rab geranylgeranyltransferase, beta subunit	
Isotype: IgG	Calculated MW: 331 aa, 37 kDa	
Immunogen Catalog Number: AG2376	Observed MW: 35 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11801-1-AP (RABGGTB antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.