

For Research Use Only

# RGS18 Polyclonal antibody, PBS Only

Catalog Number: 11866-1-PBS



## Basic Information

### Catalog Number:

11866-1-PBS

### Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG2445

### GenBank Accession Number:

BC020632

### GeneID (NCBI):

64407

### UNIPROT ID:

Q9NS28

### Full Name:

regulator of G-protein signaling 18

### Calculated MW:

235 aa, 28 kDa

### Observed MW:

28 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

WB, Indirect ELISA

### Species Specificity:

human

## Storage

### Storage:

Store at -80°C.

### Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human peripheral blood leukocyte sample were subjected to SDS PAGE followed by western blot with 11866-1-AP (RGS18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11866-1-PBS in a different storage buffer formulation.