

For Research Use Only

NUMB Polyclonal antibody

Catalog Number: 11876-1-AP



Basic Information

Catalog Number:

11876-1-AP

Size:

150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2465

GenBank Accession Number:

BC020788

GeneID (NCBI):

8650

UNIPROT ID:

P49757

Full Name:

numb homolog (Drosophila)

Calculated MW:

66 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

NUMB plays a role in the process of neurogenesis. It is required throughout embryonic neurogenesis to maintain neural progenitor cells, also called radial glial cells (RGCs), by allowing their daughter cells to choose progenitor over neuronal cell fate. NUMB may also mediate local repair of brain ventricular wall damage. NUMB may play a role in further understanding the NUMB function in development and cancer.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

