

For Research Use Only

OGFRL1 Polyclonal antibody, PBS Only

Catalog Number: 11902-1-PBS



Basic Information

Catalog Number:

11902-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2518

GenBank Accession Number:

BC018688

GeneID (NCBI):

79627

UNIPROT ID:

Q5TC84

Full Name:

opioid growth factor receptor-like 1

Calculated MW:

450 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11902-1-AP (OGFRL1 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11902-1-PBS in a different storage buffer formulation.