

For Research Use Only

# SNF8 Polyclonal antibody

Catalog Number: 11980-1-AP

Featured Product

1 Publications



## Basic Information

### Catalog Number:

11980-1-AP

### Size:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG2606

### GenBank Accession Number:

BC008976

### GeneID (NCBI):

11267

### UNIPROT ID:

Q96H20

### Full Name:

SNF8, ESCRT-II complex subunit, homolog (S. cerevisiae)

### Calculated MW:

258 aa, 29 kDa

### Observed MW:

30 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Cited Applications:

WB

### Species Specificity:

human, mouse, rat

### Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Kattareeya Kumthip	29084253	PLoS Pathog	WB

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

