

For Research Use Only

XAGE1A Polyclonal antibody

Catalog Number: 12061-1-AP



Basic Information

Catalog Number:

12061-1-AP

Size:

150ul , Concentration: 147 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2696

GenBank Accession Number:

BC009538

GeneID (NCBI):

653219

UNIPROT ID:

Q9HD64

Full Name:

X antigen family, member 1A

Calculated MW:

160 aa, 18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

XAGE1A, also named as XAGE1 and GAGED2, is a cancer testis antigen with an expression pattern that is limited to germ cells of the testis and a variety of neoplastic tissues. XAGE1A is translated as two proteins, a 9 kD protein, termed p9, and a 16.3 kD protein, termed p16. It may be useful as a target for many cancer immunotherapies.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

