

For Research Use Only

APIP Polyclonal antibody, PBS Only

Catalog Number:12089-2-PBS



Basic Information

Catalog Number:

12089-2-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2729

GenBank Accession Number:

BC009077

GeneID (NCBI):

51074

UNIPROT ID:

Q96GX9

Full Name:

APAF1 interacting protein

Calculated MW:

242 aa, 27 kDa

Observed MW:

27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12089-2-AP (AIP antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12089-2-PBS in a different storage buffer formulation.