

For Research Use Only

HYLS1 Polyclonal antibody

Catalog Number: 12141-1-AP



Basic Information

Catalog Number:

12141-1-AP

Size:

150ul , Concentration: 180 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2789

GenBank Accession Number:

BC015047

GeneID (NCBI):

219844

UNIPROT ID:

Q96M11

Full Name:

hydrolethalus syndrome 1

Calculated MW:

34 aa, 4 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse

Background Information

HYLS1, also named as HLS, is widely expressed in fetal tissue. Mutation of HYLS1 will cause the hydrolethalus syndrome (HLS) which is an autosomal recessive lethal malformation syndrome characterized by multiple developmental defects of fetus.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

