

For Research Use Only

HYLS1 Polyclonal antibody

Catalog Number: 12141-1-AP



Basic Information

Catalog Number: 12141-1-AP	GenBank Accession Number: BC015047	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 180 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 219844	
Source: Rabbit	Full Name: hydrolethalus syndrome 1	
Isotype: IgG	Calculated MW: 34 aa, 4 kDa	
Immunogen Catalog Number: AG2789		

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human, mouse

Background Information

HYLS1, also named as HLS, is widely expressed in fetal tissue. Mutation of HYLS1 will cause the hydrolethalus syndrome (HLS) which is an autosomal recessive lethal malformation syndrome characterized by multiple developmental defects of fetus.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data