

For Research Use Only

ATL1 Polyclonal antibody

Catalog Number: 12149-1-AP

2 Publications



Basic Information

Catalog Number:

12149-1-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2797

GenBank Accession Number:

BC010708

GeneID (NCBI):

51062

UNIPROT ID:

Q8WXF7

Full Name:

atlastin GTPase 1

Calculated MW:

64 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, ColP

Species Specificity:

human, mouse, rat

Cited Species:

human

Notable Publications

| Author | Pubmed ID | Journal | Application |
|------------------|-----------|---------------|-------------|
| Yupeng Zhao | 32916628 | iScience | WB, ColP |
| Eliska Zlamalova | 38851544 | Neurobiol Dis | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data