

For Research Use Only

# GPS2 Polyclonal antibody

Catalog Number: 12177-1-AP



---

## Basic Information

<b>Catalog Number:</b> 12177-1-AP	<b>GenBank Accession Number:</b> BC013652	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 493 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 2874	
<b>Source:</b> Rabbit	<b>Full Name:</b> G protein pathway suppressor 2	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 327 aa, 37 kDa	
<b>Immunogen Catalog Number:</b> AG2821	<b>Observed MW:</b> 37 kDa	

---

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

---

## Background Information

---

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data