

For Research Use Only

BFL-1/A1 Polyclonal antibody

Catalog Number: 12223-1-AP **2 Publications**



Basic Information

Catalog Number: 12223-1-AP	GenBank Accession Number: BC016281	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 273 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 597	
Source: Rabbit	Full Name: BCL2-related protein A1	
Isotype: IgG	Calculated MW: 175 aa, 20 kDa	
Immunogen Catalog Number: AG2844		

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yan-Wei Hu	26201458	Apoptosis	WB
Ming-Da Wang	26698170	Hepatology	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data