

For Research Use Only

RAB7B Polyclonal antibody, PBS Only

Catalog Number:12283-1-PBS



Basic Information

Catalog Number:

12283-1-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2936

GenBank Accession Number:

BC017092

GeneID (NCBI):

338382

UNIPROT ID:

Q96AH8

Full Name:

RAB7B, member RAS oncogene family

Calculated MW:

199 aa, 23 kDa

Observed MW:

27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12283-1-AP (RAB7B antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12283-1-PBS in a different storage buffer formulation.