

For Research Use Only

MTMR8 Polyclonal antibody, PBS Only

Catalog Number:12299-1-PBS



Basic Information

Catalog Number:

12299-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2944

GenBank Accession Number:

BC012399

GeneID (NCBI):

55613

UNIPROT ID:

Q96EFO

Full Name:

myotubularin related protein 8

Calculated MW:

704 aa, 79 kDa

Observed MW:

55-58 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 12299-1-AP (MTMR8 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12299-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of A375 cells using 12299-1-AP (MTMR8 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 12299-1-PBS in a different storage buffer formulation.