

For Research Use Only

EGFL7 Polyclonal antibody

Catalog Number: 12317-1-AP



Basic Information

Catalog Number:
12317-1-AP

Size:
150ul , Concentration: 320 µg/ml by
Bradford method using BSA as the
standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG2970

GenBank Accession Number:
BC012377

GeneID (NCBI):
51162

UNIPROT ID:
Q9UHF1

Full Name:
EGF-like-domain, multiple 7

Calculated MW:
273 aa, 30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data