

For Research Use Only

# PLEKHA4 Polyclonal antibody

Catalog Number: 12346-1-AP



## Basic Information

<b>Catalog Number:</b> 12346-1-AP	<b>GenBank Accession Number:</b> BC024157	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 173 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 57664	
<b>Source:</b> Rabbit	<b>Full Name:</b> pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 4	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 779 aa, 85 kDa	
<b>Immunogen Catalog Number:</b> AG3000		

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data