

For Research Use Only

QTRT1 Polyclonal antibody, PBS Only

Catalog Number: 12349-1-PBS

Featured Product



Basic Information

Catalog Number:

12349-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3010

GenBank Accession Number:

BC015350

GeneID (NCBI):

81890

UNIPROT ID:

Q9BXRO

Full Name:

queuine tRNA-ribosyltransferase 1

Calculated MW:

403 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of QTRT1 antibody (12349-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-QTRT1 transfected HeLa cells. This data was developed using the same antibody clone with 12349-1-PBS in a different storage buffer formulation.