

For Research Use Only

# IL-17RA Polyclonal antibody

Catalog Number: 12402-1-AP



## Basic Information

**Catalog Number:**

12402-1-AP

**Size:**

150ul , Concentration: 173 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG2955

**GenBank Accession Number:**

BC011624

**GeneID (NCBI):**

23765

**UNIPROT ID:**

Q96F46

**Full Name:**

interleukin 17 receptor A

**Calculated MW:**

866 aa, 96 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human, mouse, rat

## Background Information

IL17RA, also named as CD217 and IL17R, is a receptor for IL17A and IL17F. It is a ubiquitously expressed receptor that is essential for IL-17 biologic activity.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

