

For Research Use Only

DPM1 Polyclonal antibody, PBS Only

Catalog Number: 12403-2-PBS

Featured Product



Basic Information

Catalog Number:

12403-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1128

GenBank Accession Number:

BC007073

GeneID (NCBI):

8813

UNIPROT ID:

O60762

Full Name:

dolichyl-phosphate
mannosyltransferase polypeptide 1,
catalytic subunit

Calculated MW:

260 aa, 30 kDa

Observed MW:

27-32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 12403-2-AP (DPM1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12403-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 12403-2-AP (DPM1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12403-2-PBS in a different storage buffer formulation.