

For Research Use Only

ACRBP Polyclonal antibody, PBS Only

Catalog Number:12487-1-PBS



Basic Information

Catalog Number:
12487-1-PBS

Size:
100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG3180

GenBank Accession Number:
BC033010

GeneID (NCBI):
84519

UNIPROT ID:
Q8NEB7

Full Name:
acrosin binding protein

Calculated MW:
543 aa, 61 kDa

Observed MW:
60 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 12487-1-AP (ACRBP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12487-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 12487-1-AP (ACRBP antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12487-1-PBS in a different storage buffer formulation.