

For Research Use Only

ID2 Polyclonal antibody, PBS Only

Catalog Number: 12502-1-PBS



Basic Information

Catalog Number:

12502-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3179

GenBank Accession Number:

BC030639

GeneID (NCBI):

3398

UNIPROT ID:

Q02363

Full Name:

inhibitor of DNA binding 2, dominant negative helix-loop-helix protein

Calculated MW:

134 aa, 15 kDa

Observed MW:

15 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12502-1-AP (ID2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12502-1-PBS in a different storage buffer formulation.