

For Research Use Only

RAPGEF5 Polyclonal antibody, PBS Only

Catalog Number:12556-1-PBS



Basic Information

Catalog Number:

12556-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3229

GenBank Accession Number:

BC039203

GeneID (NCBI):

9771

UNIPROT ID:

Q92565

Full Name:

Rap guanine nucleotide exchange factor (GEF) 5

Calculated MW:

68 kDa

Observed MW:

68-75 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12556-1-AP (RAPGEF5 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12556-1-PBS in a different storage buffer formulation.



IP result of anti-RAPGEF5 (IP:12556-1-AP, 4ug; Detection:12556-1-AP 1:500) with K-562 cells lysate 1760 ug. This data was developed using the same antibody clone with 12556-1-PBS in a different storage buffer formulation.