

For Research Use Only

PDE4DIP Polyclonal antibody

Catalog Number: 12637-1-AP



Basic Information

Catalog Number:

12637-1-AP

Size:

150ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3332

GenBank Accession Number:

BC025406

GeneID (NCBI):

9659

UNIPROT ID:

Q5VU43

Full Name:

phosphodiesterase 4D interacting protein

Calculated MW:

310 aa, 36 kDa

Observed MW:

250-270 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HEK-293 cells, PC-3 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12637-1-AP (PDE4DIP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.