

For Research Use Only

# PDE4DIP Polyclonal antibody, PBS Only

Catalog Number:12637-1-PBS



## Basic Information

**Catalog Number:**  
12637-1-PBS

**Size:**  
100ug , Concentration: 1 mg/ml by  
Nanodrop;

**Source:**  
Rabbit

**Isotype:**  
IgG

**Immunogen Catalog Number:**  
AG3332

**GenBank Accession Number:**  
BC025406

**GeneID (NCBI):**  
9659

**UNIPROT ID:**  
Q5VU43

**Full Name:**  
phosphodiesterase 4D interacting  
protein

**Calculated MW:**  
310 aa, 36 kDa

**Observed MW:**  
250-270 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
human

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12637-1-AP (PDE4DIP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12637-1-PBS in a different storage buffer formulation.