

For Research Use Only

IQSEC1 Polyclonal antibody, PBS Only

Catalog Number: 12743-1-PBS



Basic Information

Catalog Number:

12743-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3681

GenBank Accession Number:

BC010267

GeneID (NCBI):

9922

UNIPROT ID:

Q6DN90

Full Name:

IQ motif and Sec7 domain 1

Calculated MW:

814 aa, 92 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data