

For Research Use Only

DNAJA4 Polyclonal antibody

Catalog Number: 12806-1-AP



Basic Information

Catalog Number:

12806-1-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3558

GenBank Accession Number:

BC021720

GeneID (NCBI):

55466

UNIPROT ID:

Q8WW22

Full Name:

DnaJ (Hsp40) homolog, subfamily A, member 4

Calculated MW:

397 aa, 45 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 12806-1-AP (DNAJA4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.