

For Research Use Only

# MAPK6 Polyclonal antibody

Catalog Number: 12839-1-AP

Featured Product

3 Publications



## Basic Information

<b>Catalog Number:</b> 12839-1-AP	<b>GenBank Accession Number:</b> BC035492	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 400 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 5597	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>Full Name:</b> mitogen-activated protein kinase 6	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 83 kDa	
<b>Immunogen Catalog Number:</b> AG3577	<b>Observed MW:</b> 100-105 kDa	

## Applications

<b>Tested Applications:</b> WB,ELISA	<b>Positive Controls:</b> WB : A431 cells,
<b>Cited Applications:</b> IF, IHC, WB	
<b>Species Specificity:</b> human, mouse, rat	
<b>Cited Species:</b> human, rat	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Guanqun Huang	29764463	J Exp Clin Cancer Res	WB,IHC,IF
Chuanxian Hu	29456679	Exp Ther Med	WB
Juanjuan Tan	28429763	Sci Rep	WB,IF

## Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 12839-1-AP (MAPK6 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.