

For Research Use Only

PAMCI Polyclonal antibody

Catalog Number: 12861-1-AP



Basic Information

Catalog Number: 12861-1-AP	GenBank Accession Number: BC031589	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop;	GeneID (NCBI): 9182	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: Ras association (RalGDS/AF-6) domain family (N-terminal) member 9	
Isotype: IgG	Calculated MW: 435 aa, 50 kDa	
Immunogen Catalog Number: AG3976	Observed MW: 50-55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, human liver tissue
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 12861-1-AP (PAMCI antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.