

For Research Use Only

ALOX5AP Polyclonal antibody

Catalog Number: 12864-1-AP



Basic Information

Catalog Number:

12864-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3144

GenBank Accession Number:

BC018538

GeneID (NCBI):

241

UNIPROT ID:

P20292

Full Name:

arachidonate 5-lipoxygenase-activating protein

Calculated MW:

18 kDa

Observed MW:

15-17 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : U-937 cells, human saliva tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12864-1-AP (ALOX5AP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.