

For Research Use Only

SECTM1 Polyclonal antibody

Catalog Number: 12865-1-AP



Basic Information

Catalog Number: 12865-1-AP	GenBank Accession Number: BC017716	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 6398	
Source: Rabbit	Full Name: secreted and transmembrane 1	
Isotype: IgG	Calculated MW: 248 aa, 27 kDa	
Immunogen Catalog Number: AG3190		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

SECTM1, also named as protein K12, is identified by its location just upstream of the CD7 locus. It's a novel molecule with potential importance in hematopoietic and/or immune system processes. The calculated MW of SECTM1 is 27 kDa. It can be cleaved into a 20 kDa soluble form.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data