

For Research Use Only

RICTOR Polyclonal antibody

Catalog Number: 12901-1-AP

1 Publications



Basic Information

Catalog Number:

12901-1-AP

Size:

150ul , Concentration: 180 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3598

GenBank Accession Number:

BC029608

GeneID (NCBI):

253260

UNIPROT ID:

Q6R327

Full Name:

rapamycin-insensitive companion of mTOR

Calculated MW:

192 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Hui Wang	25605807	Diabetes	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 12901-1-AP (RICTOR antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.