

For Research Use Only

# UPP2 Polyclonal antibody, PBS Only

Catalog Number: 12964-1-PBS



## Basic Information

**Catalog Number:**

12964-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by  
Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG3619

**GenBank Accession Number:**

BC033529

**GeneID (NCBI):**

151531

**UNIPROT ID:**

O95045

**Full Name:**

uridine phosphorylase 2

**Calculated MW:**

317 aa, 36 kDa

**Observed MW:**

36 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

IP, Indirect ELISA

**Species Specificity:**

human, mouse

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP result of anti-UPP2 (IP:12964-1-AP, 4µg; Detection:12964-1-AP 1:500) with mouse kidney tissue lysate 3400 µg. This data was developed using the same antibody clone with 12964-1-PBS in a different storage buffer formulation.