

For Research Use Only

PASD1 Polyclonal antibody

Catalog Number: 12976-1-AP



Basic Information

Catalog Number: 12976-1-AP	GenBank Accession Number: BC040301	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 750 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 139135	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	Full Name: PAS domain containing 1	
Isotype: IgG	Calculated MW: 773 aa, 87 kDa	
Immunogen Catalog Number: AG3587	Observed MW: 72-87, 150 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : NCI-H1299 cells, U2OS cells
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12976-1-AP (PASD1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.