

For Research Use Only

MSL3 Polyclonal antibody

Catalog Number: 12977-1-AP



Basic Information

Catalog Number:

12977-1-AP

Size:

150ul , Concentration: 1100 ug/ml by Nanodrop and 600 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3592

GenBank Accession Number:

BC031210

GeneID (NCBI):

10943

UNIPROT ID:

Q8N5Y2

Full Name:

male-specific lethal 3 homolog (Drosophila)

Calculated MW:

521 aa, 60 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse skeletal muscle tissue, human stomach tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 12977-1-AP (MSL3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.