

For Research Use Only

RHOXF2 Polyclonal antibody

Catalog Number: 13088-1-AP

1 Publications



Basic Information

Catalog Number:

13088-1-AP

Size:

150ul , Concentration: 1000 ug/ml by Nanodrop and 400 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3602

GenBank Accession Number:

BC021719

GeneID (NCBI):

84528

UNIPROT ID:

Q9BQY4

Full Name:

Rhox homeobox family, member 2

Calculated MW:

32 kDa

Observed MW:

35-38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : K-562 cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Jennifer Borgmann	28171660	Hum Mol Genet	WB,IHC,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 13088-1-AP (RHOXF2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.