

For Research Use Only

CHST9 Polyclonal antibody, PBS Only

Catalog Number: 13135-1-PBS



Basic Information

Catalog Number:

13135-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3874

GenBank Accession Number:

BC025764

GeneID (NCBI):

83539

UNIPROT ID:

Q7L155

Full Name:

carbohydrate (N-acetylgalactosamine 4-O) sulfotransferase 9

Calculated MW:

438 aa, 52 kDa

Observed MW:

52 kDa, 42 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 13135-1-AP (CHST9 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13135-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human intrahepatic cholangiocarcinoma tissue slide using 13135-1-AP (CHST9 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 13135-1-PBS in a different storage buffer formulation.