

For Research Use Only

ZC3HC1 Polyclonal antibody, PBS Only

Catalog Number: 13148-1-PBS



Basic Information

Catalog Number:

13148-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3744

GenBank Accession Number:

BC028917

GeneID (NCBI):

51530

UNIPROT ID:

Q86WB0

Full Name:

zinc finger, C3HC-type containing 1

Calculated MW:

502 aa, 55 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cell lysates were subjected to SDS PAGE followed by western blot with 13148-1-AP (ZC3HC1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13148-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 13148-1-AP (ZC3HC1 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 13148-1-PBS in a different storage buffer formulation.